

---

## **Hi-Torque<sup>®</sup> Guide Wires**

### **CAUTION**

Federal (USA) law restricts this device to sale by or on the order of a physician.

This device should be used only by physicians trained in angiography and percutaneous transluminal coronary angioplasty (PTCA), and/or percutaneous transluminal angioplasty (PTA).

CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. OBSERVE ALL WARNINGS AND PRECAUTIONS NOTED THROUGHOUT THESE INSTRUCTIONS. FAILURE TO DO SO MAY RESULT IN COMPLICATIONS.

Refer to the instructions supplied with any interventional devices to be used in conjunction with the Hi-Torque Guide Wire for their intended uses, contraindications, and potential complications.

### **DESCRIPTION**

The Hi-Torque Guide Wire is a steerable guide wire available in several lengths and diameters. The distal tip is shapeable or, as an option, a preshaped “J” tip is available. Refer to the product label for product specifications (e.g., wire length, diameter, and length of tip radiopacity).

Hi-Torque Extendable Guide Wires: Some Hi-Torque Guide Wires have a modified proximal end that permits the attachment of the DOC<sup>®</sup> Guide Wire Extension. Refer to the product label for Guide Wire Extension system compatibility. Joining the Guide Wire Extension to the guide wire facilitates the exchange of one interventional device for another, while maintaining guide wire position in the anatomy. After the interventional device exchange has been completed, the extension can be detached and the guide wire can be used in its original capacity.

.014 Hi-Torque Guide Wires with Proximal Markers: Brachial and femoral markers located on the proximal segment of the 0.014” (0.36 mm) guide wire aid in gauging guide wire position relative to the guiding catheter tip when using bare wire technique. They are compatible with guiding catheters that are at least 90 cm (brachial) or 100 cm (femoral) long.

Hi-Torque Guide Wires with Distal Radiopaque Markers: Some Hi-Torque Guide Wires have distal radiopaque markers. Refer to the product drawing on the label to determine the presence and location of the markers. Markers are represented by the following designation: ■. A series of spaced markers are provided proximal to the coil as references in determining lesion length.

Hi-Torque Guide Wires with Hydrophilic Coating: Refer to the product label for presence of a hydrophilic coating. When wet, a hydrophilic coating increases the lubricity of the guide wire surface.

---

## HOW SUPPLIED

**Sterile.** Sterilized with ethylene oxide gas or electron-beam radiation. Refer to the product label for the sterile method. Non-pyrogenic. Do not use if the package is open or damaged.

**Contents.** One (1) Guide Wire

**Storage.** Store in a dry, dark, cool place.

## INTENDED USE

All Hi-Torque Guide Wires are intended to facilitate the placement of balloon dilatation catheters during percutaneous transluminal coronary angioplasty (PTCA) and percutaneous transluminal angioplasty (PTA).

## LABEL DESIGNATIONS

Designations	Meaning
STENT	Intravascular stents
ATHERECTOMY	Intravascular directional atherectomy devices
IVUS	Intravascular ultrasound devices

## INDICATIONS

Refer to the device label for any additional product specific indications that may apply.

## CONTRAINDICATIONS

Hi-Torque Guide Wires are not intended for use in the cerebral vasculature.

Refer to the device label for any additional product specific contraindications that may apply.

## WARNINGS

This device is designed and intended for ONE TIME USE ONLY. DO NOT RESTERILIZE AND / OR REUSE.

Observe guide wire movement in the vessels. Before a guide wire is moved or torqued, the tip movement should be examined under fluoroscopy. Do not torque a guide wire without observing corresponding movement of the tip; otherwise, vessel trauma may occur. In addition, during catheter manipulations, ensure that the distal guide wire tip is visible.

---

Torquing a guide wire against resistance may cause guide wire damage and/or guide wire tip separation. Always advance or withdraw the guide wire slowly. Never push, auger, withdraw or torque a guide wire, which meets resistance. Resistance may be felt and / or observed under fluoroscopy by noting any buckling of the guide wire tip. If guide wire tip prolapse is observed or used for positioning, do not allow the tip to remain in a prolapsed condition; otherwise, damage to the guide wire may occur. Determine the cause of resistance under fluoroscopy and take any necessary remedial action.

If the guide wire tip becomes entrapped within the vasculature, **DO NOT TORQUE THE GUIDE WIRE.**

Maintain continuous flush while removing and reinserting the guide wire to prevent air from entering the catheter system. Perform exchanges slowly to prevent air entry and / or trauma.

When reintroducing the guide wire, confirm that the interventional device tip is free within the vessel lumen and not against the vessel wall. Failure to do so may result in vessel trauma upon guide wire exit of the device. Use the radiopaque marker of the interventional device to confirm position.

## **PRECAUTIONS**

Guide wires are delicate instruments and should be handled carefully. Prior to use and when possible during the procedure, inspect the guide wire carefully for bends, kinks, or other damage. Do not use damaged guide wires. Using a damaged guide wire may result in vessel damage and/or inaccurate torque response.

Confirm the compatibility of the guide wire diameter with the interventional device before actual use.

Free movement of the guide wire within the interventional device is an important feature of a steerable guide wire system because it gives the user valuable tactile information. Test the system for any resistance prior to use. Adjust or replace the hemostatic valve with an adjustable valve if it is found to inhibit guide wire movement.

Never attach the torque device to the modified portion of the proximal end of the extendable guide wire; otherwise, guide wire damage may occur, preventing the ability to attach the DOC Guide Wire Extension.

Hi-Torque Guide Wires with Hydrophilic Coating: Avoid abrasion of the hydrophilic coating.

Do not withdraw or manipulate the hydrophilic-coated wire in a metal cannula or sharp-edged object.

## **PREPARATION FOR USE**

Contraindications, warnings and intended use of interventional devices compatible with Hi-Torque Guide Wires are described in the instructions supplied with the respective device.

Prior to the interventional procedure, examine carefully all equipment to be used, including the interventional device, for defects. Do not use any defective equipment.

- 
1. Prepare the interventional device according to the manufacturer's instructions. Be sure to flush the guide wire lumen before introducing the guide wire.
  2. Remove the guide wire from the dispenser by pushing the exposed section of the guide wire into the dispenser until the guide wire tip and a portion of the core exit the end of the hoop. Then grasp the core of the wire to remove it totally from the dispenser. Avoid damaging the fragile guide wire tip. Do not grasp the tip of the wire while removing it from the dispenser.
  3. If indicated, the guide wire tip may be carefully shaped using standard tip shaping practices. Do not use a shaping instrument with a sharp edge.

### **Hi-Torque Guide Wires with Hydrophilic Coating**

1. Before removing the guide wire from the wire dispenser, inject normal saline into the hub end of the dispenser to thoroughly wet the complete surface of the guide wire.
2. Carefully remove the guide wire from the dispenser, as suggested above in Preparation For Use, Step 2. If the guide wire cannot be removed easily from the dispenser, inject more normal saline and attempt to remove the guide wire again.
3. Do not reinsert the guide wire into the dispenser once it has been removed.
4. If the surface of the hydrophilic-coated wire becomes dry, wetting the surface with normal saline will renew the hydrophilic effect. Be sure to thoroughly rewet the guide wire before reintroduction into an interventional device.
5. After the guide wire is withdrawn from the body, it should be wiped clean with saline-soaked gauze and kept wet.

## **DIRECTIONS FOR USE**

### **Over-The-Wire Type Systems (Preload Technique)**

1. Carefully insert the guide wire through the guide wire lumen hub of the interventional device.
2. Advance the guide wire until its tip is just proximal to the interventional device tip.
3. If using a guiding catheter, engage the guiding catheter and insert the interventional device / guide wire assembly through the hemostatic valve. Advance the system through the guiding catheter until it is just proximal to the tip of the guiding catheter.
4. Tighten the hemostatic valve to create a seal around the interventional device. Ensure intentional guide wire movement is still permitted.
5. Attach the torque device to the guide wire, if desired.
6. Under fluoroscopy, advance the guide wire out of the interventional device while securing the interventional device in place. Use the torque device to steer the guide wire across the lesion.

- 
7. Secure the guide wire in place while tracking the interventional device over it and into the lesion.
  8. If a different tip configuration or guide wire is indicated, carefully remove the guide wire while observing guide wire movement under fluoroscopy.
  9. Reshape the guide wire tip according to standard practice or prepare the next guide wire to be used.
  10. Reinsert the guide wire following Steps 1 through 7 of this section.

### **Rail Type Systems (Bare Wire Technique)**

1. Engage the guiding catheter and then insert a guide wire introducer through the hemostatic valve attached to the guiding catheter.
2. Carefully insert the distal tip of the guide wire through the introducer and into the guiding catheter.
3. Hi-Torque Guide Wires with Hydrophilic Coating: If a metal guide wire introducer was used, be sure to remove it before withdrawing or further manipulating the wire.
4. .014 Hi-Torque Guide Wires with Proximal Markers: Advance the guide wire to the appropriate proximal marker. When the proximal marker is aligned with the knurled knob of the hemostatic valve, the guide wire tip is just proximal to the guiding catheter tip.

**NOTE.** Use the most distal marker as a distance gauge when using a 90 cm brachial guiding catheter, and the most proximal marker as the distance gauge when using a 100 cm femoral guiding catheter.

5. Attach the torque device.
6. Under fluoroscopy, advance the guide wire out of the guiding catheter and into the selected vessel. Use the torque device to steer the guide wire across the lesion.
7. If a different tip configuration or guide wire is indicated, the guide wire may be removed as follows:
  - a. Open the hemostatic valve and the flush line on the coronary manifold. Slowly withdraw the guide wire while observing guide wire movement under fluoroscopy.
  - b. Close the hemostatic valve and coronary manifold flush line.
8. Reshape the guide wire tip according to standard practice or prepare the next guide wire.
9. Reinsert the guide wire following Steps 2 through 6 of this section.

- 
10. Remove the torque device and the guide wire introducer from the guide wire.
  11. Secure the guide wire while tracking the interventional device over it and into the lesion.

## **INSTRUCTIONS FOR INTERVENTIONAL DEVICE EXCHANGE PROCEDURE**

### **Over-The-Wire Type Systems (Preload Technique)**

1. If using a Hi-Torque Exchange Guide Wire, proceed to Step 3.
2. If using a Hi-Torque Extendable Guide Wire, extend the guide wire using the DOC Guide Wire Extension. Refer to the instructions included with the extension for specific information on the use of the Hi-Torque Guide Wire as an exchange length guide wire.
3. Maintain guide wire position while withdrawing the interventional device over the exchange length guide wire.
4. Prepare the other interventional device per manufacturer's instructions.
5. Load the device onto the guide wire and advance it over the Hi-Torque Exchange Guide Wire and across the lesion.
6. Proceed according to standard medical practices.

### **Rail Type Systems (Bare Wire Technique)**

1. Maintain guide wire position while withdrawing the device over the guide wire.
2. Prepare the next interventional device per manufacturer's instructions.
3. Load the device onto the guide wire and advance it over the Hi-Torque Guide Wire and across the lesion.
4. Proceed according to standard medical practices.

## **PATENTS**

This product and/or its use are protected by one or more of the following United States patents: RE 34,466; 4,875,489; 5,341,818; 5,411,476; 5,507,301; 6,165,292; 6,221,425; 6,287,258; 6,379,369; 6,464,650; 6,491,648; 6,602,228; 6,673,025. Additional patents pending.

Hi-Torque and DOC are registered trademarks of Abbott Laboratories in the USA and in selected other countries.

**Abbott Vascular**

Santa Clara, CA 95054-2807 USA

**CUSTOMER SERVICE**











TEL: (800) 227-9902

FAX: (800) 601-8874

Outside USA TEL: (951) 914 -4669

Outside USA FAX: (951) 914 -2531

**Graphical Symbols for Medical Device Labeling**

 Manufacturer	 Inner Diameter
<b>REF</b> Catalogue Number	 Outer Diameter
<b>F</b> French Size	 Stent Length
 Guiding Catheter	 Date of Manufacture
 Consult Instructions For Use	 Use By
 Contents (Numeral represents quantity of units inside.)	<b>LOT</b> Batch Code
 Do Not Reuse	<b>STERILE R</b> Sterilized Using Irradiation
<b>STERILE EO</b> Sterilized Using Ethylene Oxide	