


## MULTI-LINK MINI VISION™ OTW Coronary Stent System

Stock Number	UPN Number	Stent Sizes Available (mm)	Stent Length (mm)	Stent Nominal Pressure (atm)	Rated Burst Pressure (atm)
1010152-08	00802526331008	2.00	8	8	16
1010153-08	00802526331060	2.25	8	9	16
1010154-08	00802526331121	2.50	8	9	16
1010152-12	00802526331015	2.00	12	8	16
1010153-12	00802526331077	2.25	12	9	16
1010154-12	00802526331138	2.50	12	9	16
1010152-15	00802526331022	2.00	15	8	16
1010153-15	00802526331084	2.25	15	9	16
1010154-15	00802526331145	2.50	15	9	16
1010152-18	00802526331039	2.00	18	8	16
1010153-18	00802526331091	2.25	18	9	16
1010154-18	00802526331152	2.50	18	9	16
1010152-23	00802526331046	2.00	23	8	16
1010153-23	00802526331107	2.25	23	9	16
1010154-23	00802526331169	2.50	23	9	16
1010152-28	00802526331053	2.00	28	8	16
1010153-28	00802526331114	2.25	28	9	16
1010154-28	00802526331176	2.50	28	9	16

Minimum I.D. of Guiding Catheter = .056"

Maximum Guide Wire = .014"



## MULTI-LINK MINI VISION™ RX Coronary Stent System

Stock Number	UPN Number	Stent Sizes Available (mm)	Stent Length (mm)	Stent Nominal Pressure (atm)	Rated Burst Pressure (atm)
1007821-08	00802526365553	2.00	8	8	16
1007822-08	00802526365614	2.25	8	9	16
1007823-08	00802526365676	2.50	8	9	16
1007821-12	00802526365560	2.00	12	8	16
1007822-12	00802526365621	2.25	12	9	16
1007823-12	00802526365683	2.50	12	9	16
1007821-15	00802526365577	2.00	15	8	16
1007822-15	00802526365638	2.25	15	9	16
1007823-15	00802526365690	2.50	15	9	16
1007821-18	00802526365584	2.00	18	8	16
1007822-18	00802526365645	2.25	18	9	16
1007823-18	00802526365706	2.50	18	9	16
1007821-23	00802526365591	2.00	23	8	16
1007822-23	00802526365652	2.25	23	9	16
1007823-23	00802526365713	2.50	23	9	16
1007821-28	00802526365607	2.00	28	8	16
1007822-28	00802526365669	2.25	28	9	16
1007823-28	00802526365720	2.50	28	9	16

Minimum I.D. of Guiding Catheter = .056"

Maximum Guide Wire = .014"