

OMNILINK®

Biliary Stent System

.035

STOCK NUMBER		SUFFIX STENT LENGTH (mm)	STENT DIAMETER (mm)	MAXIMUM POST DILATATION DIAMETER (mm)	RATED BURST PRESSURE (atm)	MAXIMUM CRIMPED STENT OD (mm)	MINIMUM SHEATH SIZE [ID(mm)]	MINIMUM GUIDING CATHETER [ID(in)]	MAXIMUM GUIDE WIRE (in)
PREFIX CATHETER LENGTH (cm)	80								
1008176	1008161	-12, -16, -18, -28	5.0	8.0	14	2.11	6F (2.13)	8F (.088)	.035"
1008178	1008163	-12, -16, -18, -28	6.0	8.0	14	2.11	6F (2.13)	8F (.088)	
1008180	1008165	-12, -16, -18, -28	7.0	8.0	14	2.11	6F (2.13)	8F (.088)	
1008182	1008167	-18, -28	8.0	11.0	12	2.21	7F (2.34)	8F (.088)	
1008183	1008168	-18, -28	9.0	11.0	12	2.31	7F (2.34)	9F (.098)	
1008184	1008169	-18, -28	10.0	11.0	12	2.39	8F (2.67)	9F (.098)	
1008176	1008161	-38	5.0	8.0	14	2.11	6F (2.13)	8F (.088)	
1008178	1008163	-38	6.0	8.0	14	2.11	6F (2.13)	8F (.088)	
1008180	1008165	-38	7.0	8.0	14	2.16	7F (2.34)	8F (.088)	
1008182	1008167	-38	8.0	11.0	12	2.21	7F (2.34)	8F (.088)	
1008183	1008168	-38	9.0	11.0	12	2.31	7F (2.34)	9F (.098)	
1008184	1008169	-38	10.0	11.0	12	2.39	8F (2.67)	9F (.098)	
1008176	1008161	-58	5.0	8.0	14	2.11	6F (2.13)	8F (.088)	
1008178	1008163	-58	6.0	8.0	14	2.11	6F (2.13)	8F (.088)	
1008180	1008165	-58	7.0	8.0	14	2.16	7F (2.34)	8F (.088)	
1008182	1008167	-58	8.0	11.0	12	2.31	7F (2.34)	9F (.098)	
1008183	1008168	-58	9.0	11.0	12	2.31	8F (2.67)	9F (.098)	
1008184	1008169	-58	10.0	11.0	12	2.39	8F (2.67)	9F (.098)	