



# VOYAGER™ OTW

## Coronary Dilatation Catheter

Stock Number	UPN Number	Inflated Balloon Diameter (mm)	Crossing Profile (in)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Proximal-Distal Shaft Diameters (F)
<b>8 mm</b>						
1009439-08	00802526366123	2.00	.024	8	14	3.2/2.4
1009440-08	00802526366154	2.25	.024	8	14	3.2/2.5
1009441-08	00802526366185	2.50	.024	8	14	3.2/2.5
1009442-08	00802526366215	2.75	.024	8	14	3.2/2.5
1009443-08	00802526366246	3.00	.024	8	14	3.2/2.5
1009444-08	00802526366277	3.25	.026	8	14	3.2/2.6
1009445-08	00802526366307	3.50	.026	8	14	3.2/2.6
1009446-08	00802526366338	3.75	.026	8	14	3.2/2.7
1009447-08	00802526366369	4.00	.029	8	14	3.2/2.7
<b>12 mm</b>						
1009438-12	00802526366109	1.50	.024	8	14	3.2/2.4
1009439-12	00802526366130	2.00	.024	8	14	3.2/2.4
1009440-12	00802526366161	2.25	.024	8	14	3.2/2.5
1009441-12	00802526366192	2.50	.024	8	14	3.2/2.5
1009442-12	00802526366222	2.75	.024	8	14	3.2/2.5
1009443-12	00802526366253	3.00	.024	8	14	3.2/2.5
1009444-12	00802526366284	3.25	.026	8	14	3.2/2.6
1009445-12	00802526366314	3.50	.026	8	14	3.2/2.6
1009446-12	00802526366345	3.75	.026	8	14	3.2/2.7
1009447-12	00802526366376	4.00	.029	8	14	3.2/2.7
<b>15 mm</b>						
1009438-15	00802526380709	1.50	.024	8	14	3.2/2.4
1009439-15	00802526380716	2.00	.024	8	14	3.2/2.4
1009440-15	00802526380723	2.25	.024	8	14	3.2/2.5
1009441-15	00802526380730	2.50	.024	8	14	3.2/2.5
1009442-15	00802526380747	2.75	.024	8	14	3.2/2.5
1009443-15	00802526380754	3.00	.024	8	14	3.2/2.5
1009444-15	00802526380761	3.25	.026	8	14	3.2/2.6
1009445-15	00802526380778	3.50	.026	8	14	3.2/2.6
1009446-15	00802526380785	3.75	.026	8	14	3.2/2.7
1009447-15	00802526380792	4.00	.029	8	14	3.2/2.7
<b>20 mm</b>						
1009439-20	00802526380808	2.00	.024	8	14	3.2/2.4
1009440-20	00802526380815	2.25	.024	8	14	3.2/2.5
1009441-20	00802526380822	2.50	.024	8	14	3.2/2.5
1009442-20	00802526380839	2.75	.024	8	14	3.2/2.5
1009443-20	00802526380846	3.00	.024	8	14	3.2/2.5
1009444-20	00802526380853	3.25	.026	8	14	3.2/2.6
1009445-20	00802526380860	3.50	.026	8	14	3.2/2.6
1009446-20	00802526380877	3.75	.026	8	14	3.2/2.7
1009447-20	00802526380884	4.00	.029	8	14	3.2/2.7

Stock Number	UPN Number	Inflated Balloon Diameter (mm)	Crossing Profile (in)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Proximal-Distal Shaft Diameters (F)
<b>30 mm</b>						
1009439-30	00802526380891	2.00	.024	8	14	3.2/2.5
1009440-30	00802526380907	2.25	.024	8	14	3.2/2.5
1009441-30	00802526380914	2.50	.024	8	14	3.2/2.5
1009442-30	00802526380921	2.75	.024	8	14	3.2/2.5
1009443-30	00802526380938	3.00	.024	8	14	3.2/2.6
1009444-30	00802526380945	3.25	.026	8	14	3.2/2.6
1009445-30	00802526380952	3.50	.026	8	14	3.2/2.7
1009446-30	00802526380969	3.75	.026	8	14	3.2/2.7
1009447-30	00802526380976	4.00	.029	8	14	3.2/2.7
Usable Catheter Working Length = 143 cm					Maximum Guide Wire = .014"	