

RX HERCULINK® ELITE™

Biliary Stent System

.014

STOCK NUMBER		SUFFIX STENT LENGTH (mm)	STENT DIAMETER (mm)	MAXIMUM POST DILATATION DIAMETER (mm)	RATED BURST PRESSURE (atm)	SHAFT DIAMETER PROXIMAL / DISTAL (F / F)	MAXIMUM CRIMPED STENT OD (in / mm)	MINIMUM SHEATH SIZE (F)	MINIMUM GUIDING CATHETER ID (F / in)
PREFIX CATHETER LENGTH (cm)									
80	135								
1011486	1011487	-12	4.0	7.0	14	3.3F / 2.9F	.067 / 1.70	5F	6F (.067")
1011489	1011490	-12	4.5	7.0		3.3F / 2.9F			
1011492	1011493	-12	5.0	7.0		3.3F / 2.9F			
1011495	1011496	-12	5.5	7.0		3.3F / 3.3F			
1011498	1011499	-12	6.0	7.0		3.6F / 3.3F			
1011501	1011502	-12	6.5	8.0		3.6F / 3.5F			
1011486	1011487	-15, -18	4.0	7.0		3.3F / 2.9F			
1011489	1011490	-15, -18	4.5	7.0		3.3F / 2.9F			
1011492	1011493	-15, -18	5.0	7.0		3.3F / 2.9F			
1011495	1011496	-15, -18	5.5	7.0		3.3F / 3.3F			
1011498	1011499	-15, -18	6.0	7.0		3.6F / 3.3F			
1011501	1011502	-15, -18	6.5	8.0		3.6F / 3.5F			
1011504	1011505	-15, -18	7.0	8.0		3.6F / 3.5F			