



Available in graduated tip stiffnesses: 50, 150, 200

HI-TORQUE PILOT™ 50

Guide Wire with Hydrocoat hydrophilic coating

Stock Number	UPN Number	Tip Radiopacity (cm)	Tip Shape	Wire Size (in)	Wire Length (cm)	List Price (\$)
1010480-H	00802526363931	3	Straight	.014	190	150.00
1010480-HJ	00802526363948	3	'J'	.014	190	150.00
1010483-H	00802526363993	3	Straight	.014	300	165.00
1010483-HJ	00802526364006	3	'J'	.014	300	165.00

HI-TORQUE PILOT™ 150

Guide Wire with Polymer Tip

Stock Number	UPN Number	Tip Radiopacity (cm)	Tip Shape	Wire Size (in)	Wire Length (cm)	List Price (\$)
1010481-H	00802526363955	3	Straight	.014	190	150.00
1010481-HJ	00802526363962	3	'J'	.014	190	150.00
1010484-H	00802526364013	3	Straight	.014	300	165.00
1010484-HJ	00802526364020	3	'J'	.014	300	165.00

HI-TORQUE PILOT™ 200

Guide Wire with Polymer Tip

Stock Number	UPN Number	Tip Radiopacity (cm)	Tip Shape	Wire Size (in)	Wire Length (cm)	List Price (\$)
1010482-H	00802526363979	3	Straight	.014	190	150.00
1010482-HJ	00802526363986	3	'J'	.014	190	150.00
1010485-H	00802526364037	3	Straight	.014	300	165.00
1010485-HJ	00802526364044	3	'J'	.014	300	165.00