


MULTI-LINK OTW ULTRA™ Coronary Stent System

Stock Number	UPN Number	Stent Sizes Available (mm)	Stent Length (mm)	Minimal I.D. of Guiding Catheter (in)	Stent Nominal Pressure (atm)	Rated Burst Pressure (atm)	List Price (\$)
1003391-13	00802526207372	4.50	13	.066	9	14	1,670.00
1003392-13	00802526207389	5.00	13	.075	9	14	1,670.00
1003391-18	00802526207419	4.50	18	.066	9	14	1,795.00
1003392-18	00802526207426	5.00	18	.075	9	14	1,795.00
1003391-28	00802526207457	4.50	28	.066	9	14	2,320.00
1003392-28	00802526207464	5.00	28	.075	9	14	2,320.00
1003391-38	00802526207495	4.50	38	.066	9	14	2,520.00
1003392-38	00802526207501	5.00	38	.075	9	14	2,520.00

Maximum Guide Wire = .014"



MULTI-LINK RX ULTRA™ Coronary Stent System

Stock Number	UPN Number	Stent Sizes Available (mm)	Stent Length (mm)	Minimal I.D. of Guiding Catheter (in)	Stent Nominal Pressure (atm)	Rated Burst Pressure (atm)	List Price (\$)
1003381-13	00802526207716	4.50	13	.066	9	14	1,670.00
1003382-13	00802526207723	5.00	13	.075	9	14	1,670.00
1003381-18	00802526207754	4.50	18	.066	9	14	1,795.00
1003382-18	00802526207761	5.00	18	.075	9	14	1,795.00
1003381-28	00802526207792	4.50	28	.066	9	14	2,320.00
1003382-28	00802526207815	5.00	28	.075	9	14	2,320.00
1003381-38	00802526207853	4.50	38	.066	9	14	2,520.00
1003382-38	00802526207860	5.00	38	.075	9	14	2,520.00

Maximum Guide Wire = .014"