

SYSTEM COMPLIANCE

Balloon Size (mm)	INFLATION PRESSURE (atm)																
	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
4.0	3.58	3.68	3.76	3.84	3.90	3.96	4.00	4.04	4.08	4.11	4.14	4.16	4.19	4.21	4.23	4.26	4.28
4.5	3.88	4.00	4.11	4.20	4.36	4.43	4.50	4.56	4.61	4.66	4.70	4.75	4.79	4.83	4.87	4.92	4.97
5.0	4.34	4.47	4.59	4.69	4.83	4.92	5.00	5.07	5.14	5.21	5.27	5.34	5.40	5.46	5.53	5.60	5.68
5.5	4.94	5.06	5.17	5.26	5.35	5.43	5.50	5.57	5.63	5.68	5.74	5.79	5.84	5.90	5.95	6.01	6.08
6.0	5.35	5.49	5.62	5.73	5.83	5.92	6.00	6.07	6.13	6.19	6.24	6.29	6.34	6.39	6.45	6.50	6.56
6.5	5.76	5.93	6.08	6.21	6.32	6.42	6.50	6.57	6.63	6.68	6.73	6.77	6.81	6.85	6.88	6.92	6.97
7.0	6.23	6.39	6.54	6.67	6.79	6.90	7.00	7.09	7.18	7.26	7.34	7.41	7.49	7.58	7.66	7.76	7.86
							Nominal						RBP				