

SPECIFICATIONS

STOCK NUMBER		Balloon Diameter (mm)	Crossing Profile (in)	Nominal Pressure	Rated Burst Pressure (atm)	Proximal/Distal Shaft Diameter (F)	Minimum Guiding Catheter [ID (in)]	Minimum Sheath Size	Maximum Guide Wire (in)	
Prefix										Suffix
Catheter Length										Balloon Length (mm)
80 cm	135 cm									
1008190	1008189	-15, -20, -30, -40	4.0	0.036	8	14	3.3/3.3-2.9	6F (.067)	4F	0.014
1008192	1008191	-15, -20, -30, -40	4.5	0.037	8	14	3.3/3.3-2.9	6F (.067)	4F	0.014
1008194	1008193	-15, -20, -30, -40	5.0	0.038	8	14	3.3/3.3-2.9	6F (.067)	5F	0.014
1008196	1008195	-15, -20, -30, -40	5.5	0.041	8	14	3.3/3.3	6F (.067)	5F	0.014
1008198	1008197	-15, -20, -30, -40	6.0	0.043	8	14	3.6/3.3	6F (.067)	5F	0.014
1008200	1008199	-15, -20, -30, -40	6.5	0.048	8	14	3.6/3.5	6F (.067)	5F	0.014
1008202	1008201	-15, -20, -30, -40	7.0	0.049	8	14	3.6/3.5	6F (.067)	5F	0.014