



MULTI-LINK ZETA™ Coronary Stent System

System Compliance (Stent I.D.)

Stent Diameter (mm)	Inflation Pressure (atm)												
	8 [†]	9	10	11	12	13	14	15	16	17	18	19	20
2.75	2.75	2.80	2.84	2.89	2.93	2.97	3.00	3.04	3.07	3.09	3.12	3.14	3.15
3.00	3.00	3.05	3.09	3.13	3.17	3.21	3.25	3.28	3.32	3.35	3.38	3.41	3.44
3.50	3.50	3.55	3.60	3.64	3.68	3.72	3.75	3.78	3.81	3.84	3.86	3.88	3.90
4.00	4.00	4.04	4.09	4.13	4.18	4.22	4.25	4.29	4.33	4.37	4.42	4.46	4.50
	Nominal								RBP				

[†] These nominal *in vitro* device specifications do not take into account lesion resistance. These measurements indicate stent INNER diameter. The stent sizing should be confirmed angiographically. Do not exceed RBP or expand the stent beyond 3.50 mm (for 2.75, 3.00 mm sizes) and 4.50 mm (for 3.50, 4.00 mm sizes).