



System Compliance

Balloon Diameter (mm)	INFLATION PRESSURE (atm)												
	6	7	8	9	10	11	12	13	14	15	16	17	18
4.0	3.84	3.92	4.00	4.08	4.15	4.23	4.30	4.38	4.45	4.52	4.59	4.65	4.72
5.0	4.85	4.93	5.00	5.07	5.13	5.19	5.25	5.31	5.36	5.40	5.45	5.49	5.53
6.0	5.79	5.90	6.00	6.10	6.20	6.29	6.38	6.47	6.56	6.64	6.72	6.80	6.87
7.0	6.71	6.86	7.00	7.14	7.28	7.43	7.56	7.70	7.84	7.98	8.11	8.24	8.38
8.0	7.64	7.82	8.00	8.18	8.36	8.54	8.72	8.90	9.08	9.26	9.43	9.61	9.79
9.0	8.56	8.78	9.00	9.23	9.46	9.69	9.93	10.18	10.43	10.68	10.94	-	-
10.0	9.44	9.72	10.00	10.29	10.59	10.89	11.20	11.52	11.84	12.17	12.51	-	-
12.0	12.00	12.43	12.87	13.31	13.76	14.21	14.67	-	-	-	-	-	-

 Nominal  RBP