

Specifications

STOCK NUMBER		Balloon Diameter (mm)	Crossing Profile (in)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Proximal/Distal Shaft Diameter (F)	Minimum Sheath Size	Minimum Guiding Catheter [ID (in)]	Maximum Guide Wire (in)	
Prefix										Balloon Length (mm)
Catheter Length										
80 cm	135 cm									
1009324	1009352	-20	4.0	0.044	8	14	4.1/3.5	5F	6F (.068≤)	0.018
1009326	1009354	-20	5.0	0.045	8	14	4.1/3.5	5F	6F (.068≤)	0.018
1009328	1009356	-20	6.0	0.047	8	14	4.1/3.5	5F	6F (.068≤)	0.018
1009330	1009358	-20	7.0	0.047	8	14	4.6/4.1	5F	6F (.068≤)	0.018
1009332	1009360	-20	8.0	0.047	8	14	4.6/4.1	5F	7F (.078≤)	0.018
1009333	1009361	-20	9.0	0.051	8	12	4.6/4.1	6F	8F (.088≤)	0.018
1009334	1009362	-20	10.0	0.060	8	12	4.6/4.1	6F	8F (.088≤)	0.018
1009324	1009352	-30,-40	4.0	0.044	8	14	4.1/3.5	5F	6F (.068≤)	0.018
1009326	1009354	-30,-40	5.0	0.045	8	14	4.1/3.5	5F	6F (.068≤)	0.018
1009328	1009356	-30,-40	6.0	0.047	8	14	4.1/3.5	5F	6F (.068≤)	0.018
1009330	1009358	-30,-40	7.0	0.047	8	14	4.6/4.1	5F	7F (.078≤)	0.018
1009332	1009360	-30,-40	8.0	0.047	8	14	4.6/4.1	5F	7F (.078≤)	0.018
1009333	1009361	-30,-40	9.0	0.051	8	12	4.6/4.1	6F	8F (.088≤)	0.018
1009334	1009362	-30,-40	10.0	0.060	8	12	4.6/4.1	6F	8F (.088≤)	0.018
1009336	1009364	-30,-40	12.0	0.060	6	8	4.6/4.1	7F	9F (.098≤)	0.018
1009324	1009352	-60	4.0	0.044	8	14	4.1/3.5	5F	7F (.078≤)	0.018
1009326	1009354	-60	5.0	0.045	8	14	4.1/3.5	5F	7F (.078≤)	0.018
1009328	1009356	-60	6.0	0.047	8	14	4.1/3.5	5F	7F (.078≤)	0.018
1009330	1009358	-60	7.0	0.047	8	14	4.6/4.1	5F	7F (.078≤)	0.018
1009332	1009360	-60	8.0	0.047	8	14	4.6/4.1	6F	8F (.088≤)	0.018
1009333	1009361	-60	9.0	0.051	8	12	4.6/4.1	6F	8F (.088≤)	0.018
1009334	1009362	-60	10.0	0.060	8	12	4.6/4.1	7F	9F (.098≤)	0.018