

Compliance									
PRESSURE (atm)	BALLOON DIAMETER (mm)								
	4	5	6	7	8	9	10	12	14
4	3.76	4.64	5.66	6.57	7.56	8.42	9.15	11.42	13.36
5	3.83	4.73	5.75	6.69	7.67	8.57	9.37	11.71	13.68
6	3.88	4.82	5.83	6.79	7.78	8.72	9.58	12.0	14.0
7	3.94	4.91	5.92	6.9	7.89	8.86	9.79	12.27	14.33
8	4.0	5.0	6.0	7.0	8.0	9.0	10.0	12.54	14.65
9	4.05	5.09	6.08	7.1	8.11	9.14	10.2	12.79	14.98
10	4.1	5.18	6.15	7.19	8.21	9.26	10.41	13.03	15.32
11	4.15	5.27	6.23	7.28	8.32	9.39	10.61	13.27	15.65
12	4.2	5.35	6.3	7.37	8.42	9.51	10.81	13.49	15.99
13	4.24	5.44	6.37	7.45	8.53	9.62	11.01	N/A	N/A
14	4.28	5.53	6.44	7.53	8.63	9.73	11.21	N/A	N/A
15	4.33	5.62	6.5	7.61	8.73	9.84	11.4	N/A	N/A
16	4.36	5.71	6.56	7.68	8.83	9.94	11.6	N/A	N/A
17	4.4	5.79	6.62	7.75	8.94	N/A	N/A	N/A	N/A
18	4.44	5.88	6.68	7.81	9.04	N/A	N/A	N/A	N/A

	<b>NOMINAL</b>
	<b>RATED BURST PRESSURE</b>