

Product Specifications

STOCK NUMBER			BALLOON DIAMETER (mm)	NOMINAL PRESSURE	RATED BURST PRESSURE (atm)	MINIMUM SHEATH SIZE	MINIMUM GUIDING CATHETER ID	MAXIMUM GUIDE WIRE (in)	
PREFIX CATHETER LENGTH	SUFFIX BALLOON LENGTH (mm)								
55 cm	80 cm	135 cm							
1009976	1010004	1010060	-20,-40,-60,-80,-100	4	8	14	5F	7F(.078")	0.035
1009977	1010005	1010061	-20,-40,-60,-80,-100	5	8	14	5F	7F(.078")	0.035
1009978	1010006	1010062	-20,-40,-60,-80	6	8	14	5F	7F(.078")	0.035
1009978	1010006	1010062	-100	6	8	14	6F	8F(.088")	0.035
1009979	1010007	1010063	-20	7	8	14	5F	7F(.078")	0.035
1009979	1010007	1010063	-40,-60,-80,-100	7	8	14	6F	8F(.088")	0.035
1009980	1010008	1010064	-20,-40,-60,-80	8	8	14	6F	8F(.088")	0.035
1009981	1010009	1010065	-20,-40,-60	9	8	12	6F	8F(.088")	0.035
1009982	1010010	1010066	-20,-40	10	8	12	6F	8F(.088")	0.035
1009982	1010010	1010066	-60	10	8	12	7F	9F(.098")	0.035
1009984	1010012	1010068	-20,-40	12	6	8	7F	9F(.098")	0.035
1009986	1010014	1010070	-20,-40	14	6	8	7F	9F(.098")	0.035