

## Specifications

STOCK NUMBER		Unconstrained Stent Diameter (mm)	Reference Bile Duct (mm)	Guide Wire Compatibility
Prefix	Suffix			
Catheter Length	Stent Length (mm)			
80 cm	120 cm			
1007548	1006516	5	3.6-4.5	.018"
1007549	1006517	6	4.3-5.4	.018"
1007550	1006518	7	5.0-6.3	.018"
1007551	1006519	8	5.7-7.3	.018"
1007552	1006520	9	6.4-8.2	.018"
1007545	1006521	10	7.1-9.1	.018"

6F Sheath / 8F Guide Catheter compatible for all stent sizes

Maximum outside diameter of system .084"