

Compliance								
PRESSURE (atm)	STENT DIAMETER (mm)							
	4.0	4.5	5.0	5.5	6.0	6.5	7.0	
8	4.00	4.50	5.00	5.50	6.00	6.5	7.00	
9	4.07	4.58	5.09	5.60	6.11	6.62	7.14	
10	4.13	4.66	5.18	5.70	6.21	6.74	7.27	
11	4.19	4.74	5.27	5.80	6.30	6.85	7.39	
12	4.25	4.81	5.35	5.88	6.39	6.96	7.50	
13	4.30	4.88	5.42	5.97	6.47	7.06	7.60	
14	4.35	4.94	5.49	6.04	6.55	7.15	7.69	

	NOMINAL
	RATED BURST PRESSURE