

## Product Specifications

STOCK NUMBER		STENT LENGTH (mm)	STENT DIAMETER (mm)	MAXIMUM POST DILATATION DIAMETER (mm)	RATED BURST PRESSURE (atm)	SHAFT DIAMETER (f)	CRIMPED STENT OD (in / mm)	MINIMUM SHEATH SIZE	MINIMUM GUIDING CATHETER ID	MAXIMUM GUIDE WIRE (in)
PREFIX	SUFFIX									
CATHETER LENGTH										
80 cm	135 cm									
1008012	1008013	-12	4.0	7.1	14	3.3/2.9	.061 / 1.55	5F	6F(.067")	.014
1008014	1008015	-12	4.5	7.1	14	3.3/2.9	.061 / 1.55	5F	6F(.067")	.014
1008016	1008017	-12	5.0	7.1	14	3.3/2.9	.063 / 1.60	5F	6F(.067")	.014
1008018	1008019	-12	5.5	7.1	14	3.3/3.3	.063 / 1.60	5F	6F(.067")	.014
1008020	1008021	-12	6.0	7.1	14	3.6/3.3	.064 / 1.63	5F	6F(.067")	.014
1008022	1008023	-12	6.5	7.1	14	3.6/3.5	.070 / 1.78	6F	7F(.078")	.014
1008012	1008013	-15, -18	4.0	8.0	14	3.3/2.9	.061 / 1.55	5F	6F(.067")	.014
1008014	1008015	-15, -18	4.5	8.0	14	3.3/2.9	.061 / 1.55	5F	6F(.067")	.014
1008016	1008017	-15, -18	5.0	8.0	14	3.3/2.9	.063 / 1.60	5F	6F(.067")	.014
1008018	1008019	-15, -18	5.5	8.0	14	3.3/3.3	.063 / 1.60	5F	6F(.067")	.014
1008020	1008021	-15, -18	6.0	8.0	14	3.6/3.3	.064 / 1.63	5F	6F(.067")	.014
1008022	1008023	-15, -18	6.5	8.0	14	3.6/3.5	.070 / 1.78	6F	7F(.078")	.014
1008024	1008025	-15, -18	7.0	8.0	14	3.6/3.5	.071/ 1.80	6F	7F(.078")	.014